

# **DEVERE TEST PAPERS**







# **Chlorine Test Papers**

These papers will indicate the level of chlorine in a solution. Use with DeVere Sanitizer Concentrate. Papers measure in graduations of 10ppm, 50ppm, 100ppm and 200ppm chlorine

### **Physical Properties (Typical)**

Appearance

White Paper Strips

#### **Directions for Use**

- 1. Dip test paper in solution for 1-2 seconds
- 2. Take test paper out of solution and compare results quickly with the color chart in the packaging

Available sizes: 100 strips

# **QD Test Papers**

These papers will indicate the level of quaternary sanitizer in a solution. Use with QDII Sanitizer and One Step Quaternary Cleaner and Disinfectant. Graduations that this product measures are: 0 ppm, 100 ppm, 200 ppm, and 400 ppm quaternary ammonium chloride.

# **Physical Properties (Typical)**

Appearance

White Paper Strips

#### **Directions for Use**

- 1. Dip test paper in solution for 1-2 seconds
- 2. Take test paper out of solution and compare the test strip color to the color chart in the packaging within 10 seconds

Available sizes: 50 strips

### **GS Test Papers**

These test papers will measure the amount of iodine in a solution. Use with . Precision Iodine Test Papers measure in graduations of 0 ppm,  $12\ 1\ / 2$  ppm,  $25\ ppm$ , and  $50\ ppm$  iodine.

### **Physical Properties (Typical)**

Appearance

White Paper Strips

#### **Directions for Use**

- 1. Dip strip into the solution being tested for 60 seconds without agitation.
- 2. Compare wet test strip color to the color chart in packaging

Available sizes: 100 strips